



# QUICK AND ACCURATE METHOD TO STUDY DEMULSIFICATION

Comparison of different demulsifiers



## INTRODUCTION

The effectiveness of a demulsifier depends on many parameters including the amount of water produced, the quality of the oil/water interface, the quality of the water produced and the separation kinetics. The usual technique used to evaluate the effectiveness of a demulsifier is the Bottle Test. In this Application Note, we propose a faster and more accurate technique to assess the quality of a demulsification and to compare the efficiency of different demulsifiers.

## METHOD

The effectiveness of 5 different demulsifiers is assessed with a quick and simple protocol: 100 ppm of demulsifier (dosed from a stock solution using an appropriate solvent) is added to the crude oil emulsion and handshaked 50 times. The crude oil emulsion evaluated was from a UK onshore source with an API of 11° and a BS&W of 60%. All testing was carried out by Croda using demulsifiers from their Kemelix™ product range.

Demulsifier 1	Polysorbate polyester
Demulsifier 2	Resin alkoxyate
Demulsifier 3	Resin alkoxyate
Demulsifier 4	Polyol alkoxyate
Demulsifier 5	Polyimine derivative

Then, Turbiscan® measurements are performed at 60°C for 60 minutes with 1 scan per minute. The following graph displays a typical Turbiscan profile obtained.




Figure 1: Evolution of the transmission signal for Demulsifier 1

## RESULTS

The graph displayed in Figure 1 shows information regarding the **amount of produced water**, the **clarity** of produced water, the **quality of water/oil interface** and the **demulsification kinetics**.

### Volume of water produced

The volume of produced water is directly linked to the thickness of the transmission peak obtained. The following graph displays the transmission peak thickness for the 5 demulsifiers evaluated:




Figure 2: Peak thickness for the 5 demulsifiers evaluated

This graph allows us to classify these demulsifiers in terms of amount of produced water (initial volume of sample is 20 ml):

Sample	Final peak thickness (mm)	Amount of water produced (ml)	Amount of water produced (% of BS&W)	WATER VOLUME
Demulsifier 5	0.6	0.3	<1	-
Demulsifier 4	6.6	3.2	27	↓ +
Demulsifier 1	21.0	10.3	86	
Demulsifier 3	22.3	10.9	91	
Demulsifier 2	23.1	11.3	94	

This table shows that Demulsifiers 4 and 5 give relatively low amounts of produced water. On the other hand, Demulsifiers 1, 2 and 3 give comparable and significant amount of water produced.

### Separation kinetics

Another important parameter in the demulsification process is the time required to complete demulsification. If we look again at Figure 2, Demulsifiers 1, 2 and 3 can be compared in terms of



separation kinetics (4 and 5 gave incomplete or no demulsification so were not considered for this analysis).




Figure 3: Peak thickness for Demulsifiers 1, 2 and 3.

Sample	Start time (min)	End time (min)
Demulsifier 3	0	15
Demulsifier 2	1	20
Demulsifier 1	24	>60

SPEED  
+  
↑  
-

These results shows that Demulsifiers 2 and 3 lead to a fast demulsification, whereas Demulsifier 1 gives a slower demulsification.

### Clarity of water produced

The clarity of water produced can also be evaluated by assessing the transmission level of the water phase at the end of the measurement.




Figure 4: Transmission level at the end of the measurement.

Considering that the transmission level of pure water is 88%, the ratio between the final transmission level and the transmission level of pure water can be calculated, giving an index from 0 to 1, 1 being the best water clarity.

Sample	Final transmission level (%)	Turbiscan water clarity index
Demulsifier 5	0	0
Demulsifier 2	20.5	0.23
Demulsifier 3	67.3	0.76
Demulsifier 4	68	0.77
Demulsifier 1	83	0.94

WATER CLARITY  
-  
↓  
+

This table shows that Demulsifier 1 gives the highest clarity index, meaning the best produced water quality.

It is interesting to note that Demulsifier 4 leads to a high clarity index, but with a very low amount of water produced. Demulsifier 3 also leads to a high clarity index, but the transmission scan is very noisy, meaning that the produced water is contaminated with small oil droplets or inorganic salts.

### Quality of the oil/water interface

Finally, Figure 4 gives additional information about the quality of the oil/water interface. This parameter can be assessed by calculating the slope of the last transmission scan at the oil/water interface. The sharper the slope, the better the interface quality (cf. Figure 4). The following table displays the value of this slope for Demulsifiers 1 to 4 (Demulsifier 5 was not included due to no demulsification taking place).


Sample	Slope
Demulsifier 2	-0.78
Demulsifier 3	-1.5
Demulsifier 4	-1.7
Demulsifier 1	-8.2

INTERFACE QUALITY  
-  
↓  
+

This table shows that Demulsifier 1 has the highest slope, meaning the best interface quality (i.e. most defined oil/water interface). Other demulsifiers lead to similar slope values, much lower than Demulsifier 1.

## SUMMARY

This Application Notes shows a quick and simple method to compare the effectiveness of different demulsifiers **in only 60 minutes\***. The following graph summarizes the results obtained.



Demulsifiers 1 and 3 display the best performance; Demulsifier 1 in terms of quality of separation (oil/water interface quality and clarity of produced water) and Demulsifier 3 in terms of separation kinetics and the volume of produced water.

\*Measurement time can be longer if sample viscosity is higher than 300 cP at 60°C.